

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



NUTRIENT BROTH/AGAR I. CECT 1

- Beef extract 5 g
- Peptone 10 g
- NaCl 5 g
- Agar powder (only for solid media) 15 g
- Distilled water 1 L

Adjust pH to 7.2.

For Bacillus strains the addition of 10 mg $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ is recommended for sporulation.

