

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



M20S MEDIUM (Malt agar with 20% sucrose) . CECT 123

- Malt extract 400 ml
- Sucrose 200 g
- Agar powder (only for solid media) 15 g
- Distilled water up to 1 L

Adjust pH to 7.0 with KOH. Sterilization 15 minutes at 121 °C.

Malt extract: Dilute brewery malt with water to 10% sugar solution (level 10 on Brix sucrose meter) and sterilize 15 minutes at 121 °C.

CECT