

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



VNIVERSITAT
DE VALÈNCIA

Spirillum sp. MEDIUM. CECT 13

- Peptone 5.0 g
- Malt extract 3.0 g
- Yeast extract 3.0 g
- Calcium lactate 10.0 g
- Distilled water 1 L

Adjust pH to 7.0 with KOH. Dispense the medium into screw-capped bottles and sterilize by autoclaving 25 minutes at 115 ° C. There is no need to filter off the precipitate formed during preparation and sterilization. For solid medium add 1% purified agar.

CECT