

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



Bifidobacterium sp. MEDIUM. CECT 152

- Casein peptone (tryptic digest) 10 g
- Meat extract 5 g
- Yeast extract 5 g
- Glucose 10 g
- K₂HPO₄ 3 g
- Tween 80 1 ml
- Tap water 1 L

Adjust pH to 6.8.

For solid media add 15 g/L of agar powder.

After sterilization, aseptically add solutions of sodium ascorbate and cysteine-HCl to a final concentration of 1% and 0.05%, respectively.

Liquid medium not freshly prepared should be heated in a steamer for 10 min before addition of reducing substances.