

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



Peptococcus sp. MEDIUM. CECT 19

- Casein peptone 10.0 g
- Beef extract 3.0 g
- Yeast extract 3.0 g
- Glucose 2.0 g
- Tween 80 1 ml
- Cysteine hydrochloride 0.5 g*
- Salt solution 40 ml
- Distilled water 960 ml

Adjust pH to 7.2. For solid media add 15 g/L of agar powder.

Salt solution:

- CaCl_2 0.2 g
- $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ 0.2 g
- K_2HPO_4 1.0 g
- KH_2PO_4 1.0 g
- NaHCO_3 10.0 g
- NaCl 2.0 g
- Distilled water 1 L

* After autoclaving add aseptically 10 ml of a 5 % (w/v) solution of cysteine hydrochloride sterilized by filtration.