

# COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València  
Catedrático Agustín Escardino, 9  
46980 Paterna (Valencia). España  
Tel.: +34 96 354 46 12  
<http://www.cect.org>



## Xylophilus ampelinus MEDIUM. CECT 190

- Yeast Extract 10.0 g
- Galactose 20.0 g
- CaCO<sub>3</sub> 20.0 g
- Agar powder (only for solid medium) 15.0 g
- Distilled water 1L

Dissolve everything but the agar, add CO<sub>3</sub>Ca (which is insoluble). Adjust pH to 7.2 when necessary. Add the agar and melt it. Dispense in tubes agitating to distribute the CO<sub>3</sub>Ca. Sterilize at 115 ° C for 20 minutes.

CECT