

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



Erythromicrobium AND Roseococcus MEDIUM. CECT 202

- NH_4Cl 0.30 g
- K_2HPO_4 0.30 g
- $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ 0.50 g
- $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ 0.05 g
- KCl 0.30 g
- Na-acetate 1.00 g
- Yeast extract 1.00 g
- Bacto peptone 1.00 g
- Trace elements solution SL-6 (see medium [196](#)) 1.00 ml
- Vitamin B_{12} 20.00 μg
- Distilled water 1 L

Adjust pH to 7.5-7.8. For solid media add 15 g/L of agar powder. After autoclaving, add vitamin B_{12} from a filter sterilized stock solution.