

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



YDC Medium. CECT 33

- Yeast extract 5 g
- Dextrose (Glucose) 20 g
- CaCO₃ 20 g
- Agar powder (only for solid media) 15 g
- Distilled water 1 L

Dissolve the ingredients. Add CaCO₃ and if necessary, adjust pH to 7.2. Dispense in tubes, shaking continuously to distribute homogeneously the CaCO₃.

CECT