

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



Caulobacter sp. MEDIUM II. CECT 45

- Yeast extract 3 g
- Peptone 10 g
- Filtered aged sea water 1 L

Adjust pH to 7.2-7.4. For solid media add 15 g of agar powder.

Alternatively, Artificial Sea Water can be used for preparing the medium.

CECT