

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



Lactobacillus homohiochii MEDIUM. CECT 51

- Tryptone 10 g
- Yeast extract 5 g
- Beef extract 2 g
- Glucose 20 g
- KH_2PO_4 0.5 g
- K_2HPO_4 0.5 g
- $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ 0.2 g
- $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ 10 mg
- $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ 7.5 mg
- Tween 80 1 ml
- Ethanol 40 ml
- Sodium acetate 20 g
- Distilled water 1 L

Dissolve all ingredients in 950 ml of distilled water. Adjust pH to 5.2. Distribute in 9.5 ml amounts in tubes. Sterilize and add aseptically 0,4 ml of ethanol and 0.1 ml of a 3 mg/ml solution of DL-mevalonic acid.