

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



SPORULATION AGAR. CECT 58

- Yeast extract 1 g
- Beef extract 1 g
- Tryptone 2 g
- FeSO₄·7H₂O Traces
- Glucose 10 g
- Agar powder 15 g
- Distilled water 1 L

Adjust pH to 7.2. For broth, eliminate agar and reduce concentration to 1/3 of given quantities.

CECT