

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



MINERAL SALTS AGAR. CECT 62

- NaNO_3 2.00 g
- $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ 0.50 g
- KCl 0.50 g
- $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ 0.01 g
- KH_2PO_4 0.14 g
- K_2HPO_4 1.20 g
- Yeast extract 0.02 g
- Agar powder 15.00 g
- Distilled water 1 L

Adjust pH to 7.2. Place a strip of sterile filter paper on cooled agar slant and inoculate on to the filter paper strip.

CECT