

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



OATMEAL AGAR. CECT 68

- Oat flakes 80 g
- Agar powder 15 g
- Distilled water 1 L

Boil the oat flakes with 700 ml of water for 20 min. Filter through cheesecloth. Add the agar and bring volume to 1 L. Autoclave at 121 °C for 20 min.

CECT