

COLECCIÓN ESPAÑOLA DE CULTIVOS TIPO (CECT)



Parc Científic Universitat de València
Catedrático Agustín Escardino, 9
46980 Paterna (Valencia). España
Tel.: +34 96 354 46 12
<http://www.cect.org>



YM-1. CECT 80

- Yeast extract 5.00 g
- Mycopeptone 10.00 g
- Adenine 0.01 g
- Uracil 0.01 g
- Succinic acid 10.00 g
- NaOH 6.00 g
- Glucose 10.00 g
- Yeast Nitrogen Base 6.70 g
- Agar powder (only for solid medium) 15.00 g
- Distilled water 1 L

Adjust pH to 6.8.

Important: Sterilize the Yeast Nitrogen Base by filtration.

Vegetable peptone No. 1 can be used in place of mycopeptone.